Dentobuff® Strip

Intended use

Dentobuff is a strip method for estimation of salivary buffering capacity which is an indicator of the risk of dental erosion and caries^{1,2}.

Simultaneously with the determination of buffering capacity, stimulated salivary secretion rate can be determined.

Principle of the method

The Dentobuff test utilises a strip supporting a pH test pad impregnated with acid. A drop of stimulated saliva is pipetted onto the pad, and the result is read after 5 min. The method differentiates stimulated salivas with low, intermediate and high buffering capacity³.

Contents

Dentobuff Strip	Cat. No. 67670
Dentobuff test strips	10 pcs
Paraffin pellets	10 pcs
Disposable pipettes	10 pcs
Instructions for use	1 pc

Warnings and precautions

- Dentobuff is for *in vitro* diagnostic use only.
- Do not use paraffin pellets for small children. The test procedure requires paraffin pellet stimulated saliva and therefore it is not suitable for testing small children.
- Do not use the product beyond the expiry date marked on the box.
- Do not use the product if you detect uneven coloration of the pH pad affixed to the strip or if an unused strip is not yellow in colour.
- Do not touch the pH test pad on the strip.

Storage

Dentobuff is stored at room temperature (18...25°C) protected from moisture. The expiry date is marked on the box.

Test procedure

A. Stimulation of saliva secretion

1. The patient should sit in an upright, relaxed position.
2. Let the patient chew a paraffin pellet for 1 min⁴. During this time, the patient may either spit out or swallow any saliva produced.
3. The patient should then continue to chew the pellet for further 5 min during which any saliva produced is collected in a suitable vessel.
4. Measure the volume of the saliva produced.
5. Calculate the saliva secretion rate (e.g. 3.5 ml/5 min = 0.7 ml/min).

B. Determination of salivary buffering capacity

1. Place a Dentobuff strip on a steady surface covered with absorbent tissue, with the pH test pad facing upwards.
2. Apply one drop of **stimulated** saliva to the pH pad, using a disposable pipette supplied in the test kit.
3. The entire pH pad should be covered with the saliva.
4. After 5 min, compare the colour of the pH pad with the colour chart below.

Colour chart

- A. The colour of the pH pad
- B. The final pH of saliva
- C. Salivary buffering capacity

- | | | |
|----------|--------------|--------|
| A. Blue | Green | Yellow |
| B. ≥ 6.0 | 4.5–5.5 | ≤ 4.0 |
| C. High | Intermediate | Low |

Note 1. The test should be performed as soon as possible after saliva stimulation.

Note 2. To avoid daily variation, samples from the same patient should be taken at the same time of day during follow-up.

Note 3. Before sampling the patient should avoid

- eating for 1–3 h
- drinking for 1–3 h
- smoking for 1–3 h
- tooth brushing for 1–3 h

Note 4. Viscous saliva may cause the colour reaction to be uneven or mixed. In such cases, buffering capacity should be evaluated according to the colour indicating the lowest value.

Expected values

In an adult population (18–65 years), values < 5.0 are encountered in 17.3 %, values in the range 5.0–5.5 in 25.3 % and values > 6.0 in 57.3 % of stimulated saliva samples⁵. Bicarbonate is the chief determinant of salivary buffering capacity. The bicarbonate concentration, in turn, is strongly dependent on salivary secretion rate⁶.

The normal salivary secretion rate in a young population is > 1 ml/min⁴.

Limitations of the procedure

Eating, drinking and medications affect salivary secretion rate and should be avoided prior to testing^{1,5,6}.

Because any factor affecting salivary secretion rate can also affect the buffering capacity of saliva, it is important to collect stimulated saliva in peaceful surroundings and to use prestimulation⁴.

Performance characteristics

Dodds M. Comparison of Dentobuff with a published method. Data on file at Orion Diagnostica. 1996⁵.

A study was performed with 76 women and 74 men, mean age 30.2 years.

The mean buffering capacities for the three Dentobuff colour groups differed significantly from one another ($P < 0.0001$). The mean buffering capacity was pH 5.51 for a blue reading; 4.41 for a green reading and 3.83 for a yellow reading. Post hoc analysis revealed that these differences were all significant.

Means table for pH				
Colour	N	Mean	SD	SE
Blue	86	5.509	.803	.087
Green	38	4.414	.591	.096
Yellow	26	3.827	.782	.153

Spearman's rank correlation coefficient was 0.71, indicating a significant ($P < 0.0001$) correlation between paired rankings according to pH data and test output.

Quality control

Quality control tests are performed on each lot of Dentobuff at the time of manufacture. Should the user wish to perform his own quality control, the following procedure is recommended:

Prepare serial dilutions of fresh Sørensen phosphate buffer 0.1 M, pH 7.0 in distilled water. Dispense 50 µl of each dilution on a separate test strip. The colour reactions should be as follows:

Dilution	Colour reaction
0.1 M	blue
0.05 M	blue
0.025 M	blue
0.0125 M	green
0.00625 M	yellowish green
0.003125 M	yellow
Distilled water	yellow

Disposal

Used Dentobuff tests can be disposed of as normal waste.

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References • Litteratur • Литература

- Ericsson Y. Clinical investigations of the salivary buffering action. *Acta Odontol Scand* 1959;17:131–65.
- Hansel Petersson G, Twetman S, Bratthall D. Evaluation of a computer program for caries risk assessment in schoolchildren. *Caries Res* 2002;36: 327–40.
- Ericson D, Bratthall D. Simplified method to estimate salivary buffering capacity. *Scand J Dent Res* 1989; 97:405–7.
- Le Bell Y, Söderling E, Karjalainen S. Effect of repeated sampling and prestimulation on saliva buffer capacity and flow rate values in children. *Scand J Dent Res* 1991;99:505–9.
- Dodds M. Comparison of Dentobuff with a published method. Data on file at Orion Diagnostica. 1996.
- Heintze U, Birkhed D, Björn H. Secretion rate and buffer effect of resting and stimulated whole saliva as a function of age and sex. *Swed Dent J* 1983; 7:227–8.

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Literatura • Referencer • Viitteet

- 1 Ericsson Y. Clinical investigations of the salivary buffering action. Acta Odontol Scand 1959;17:131–65.
- 2 Hansel Petersson G, Twetman S, Bratthall D. Evaluation of a computer program for caries risk assessment in schoolchildren. Caries Res 2002;36: 327–40.
- 3 Ericson D, Bratthall D. Simplified method to estimate salivary buffering capacity. Scand J Dent Res 1989; 97:405–7.
- 4 Le Bell Y, Söderling E, Karjalainen S. Effect of repeated sampling and prestimulation on saliva buffer capacity and flow rate values in children. Scand J Dent Res 1991;99:505–9.
- 5 Dodds M. Comparison of Dentobuff with a published method. Data on file at Orion Diagnostica. 1996.
- 6 Heintze U, Birkhed D, Björn H. Secretion rate and buffer effect of resting and stimulated whole saliva as a function of age and sex. Swed Dent J 1983; 7:227–8.

Vysvětlivky použitých symbolů • Förklaring av symboler • Symbolien selitykset

	Česky	Svenska	Suomi
	Diagnostický zdravotnický prostředek <i>in vitro</i>	För <i>in vitro</i> diagnostik	<i>In vitro</i> -diagnoosiikkaan
	Katalogové číslo	Artikelnummer	Tuotenumero
	Číslo šarže	Lot. Nr.	Eräkkoodi
	Použitelné do	Används före	Käytettävä ennen
 Store at 18...25°C	Teplotní rozmezí	Förvaringstemperatur	Säilytyslämpötila
	Výrobce	Tillverkare	Valmistaja
	Obsah	Innehåll	Sisältö
	Viz návod k použití	Se bruksanvisningen	Lue käyttöohjeet

Dentobuff® is a registered trademark of Orion Diagnostica Oy.