Intended use

Slant culture tubes for primary isolation and culture of *Mycobacterium tuberculosis* complex and mycobacteria other than *M. tuberculosis*^{1,2}.

Contents

Mycobacterium tube 1 Cat. No. 68140	Mycobacterium tube 2 Cat. No. 68141
10 Mycobacterium tube 1 Instructions for use	10 Mycobacterium tube 2 Instructions for use

Material required but not supplied

Reagents and laboratory equipment for mycobacterial culture.

Principle

Mycobacterium tube 1 and 2 contain a modification of Löwenstein-Jensen medium.

Mycobacterium tube 1 supports the growth of mycobacteria that are able to utilise glycerol as the main carbon source.

Mycobacterium tube 2 is designed for isolation and culture of pyruvate-dependent mycobacteria which may remain undetected unless a pyruvate-containing medium with decreased pH is utilised^{1,2,3,4,5,6}.

Parallel use of Mycobacterium tube 1 and 2 allows reliable detection of clinically important mycobacteria.

Typical formulation

Mycobacterium tube 1		Mycobacterium tube 2	
Mycobacterium 1 salt solution	330 ml	Mycobacterium 2 salt solution	330 ml
KH ₂ PO ₄	4.0 g/l	Pyruvic acid	12.0 g/l
MgSO ₄	0.4 g/l	KH ₂ PO ₄	4.0 g/l
Trimagnesium dicitrate	1.56 g/l	MgSO ₄	0.4 g/l
Starch	50.0 g/l	Trimagnesium dicitrate	1.56 g/l
Glycerol	20 ml	Starch	50.0 g/l
Homogenised whole egg	660 ml	Homogenised whole egg	660 ml
Malachite green oxalate 1%	10 ml	Malachite green oxalate 1%	10 ml

Warnings and precautions

- Mycobacterium tube 1 and 2 are for *IN VITRO* DIAGNOSTIC USE only.
- Do not use the product beyond the expiry date.
- Do not use the Mycobacterium tube if you detect discoloration or dehydration of the normally greenish medium, poor solidification of the medium or evidence of bacterial or fungal growth.
- Because all inoculated Mycobacterium tubes may contain potential pathogens, do not touch the medium.
- Pay special attention to the procedures and devices used in mycobacterial culture because of the biohazards associated with handling of infectious material.
- Small air bubbles sometimes present inside the culture medium do not affect the performance of the medium.

Storage

Mycobacterium tube 1 and 2 should be stored at 2...8°C, protected from light and draught. **Do not allow to freeze.** The expiry date is marked on the box.

Sample collection

It is recommended that the standard procedures described by for example the United States Centers for Disease Control (CDC) be used for primary isolation of mycobacteria⁵.

Laboratory evaluation of the sample involves

1. decontamination, digestion and concentration
2. microscopic examination for the presence of acid-fast bacilli
3. isolation of the microbe by culturing
4. identification of the microbe for susceptibility testing⁶

Reading the result

Inoculated Mycobacterium tube 1 and 2 are read within five to seven (5–7) days after inoculation and then once a week for eight (8) weeks or according to local standard procedure.

Expected values

A positive *M. tuberculosis* complex culture is always considered clinically significant.

There are in addition more than 25 species in the *Mycobacterium* genus that may cause human diseases^{2,6,7}.

Performance characteristics

According to Lue et al., 1995⁸

GROWTH % (mean recovery time/days)				
Organism	N	Mycobacterium tube 1	Mycobacterium tube 2	Löwenstein-Jensen medium
<i>M. tuberculosis</i>	60	98.6 (16.4)	95.7 (16.6)	100 (23.2)
<i>M. avium</i>	17	87.0 (19.1)	95.7 (18.4)	88.2 (28.2)
<i>M. gordonae</i>	1	100 (13.0)	100 (13.0)	100 (20.0)
<i>M. kansasii</i>	1	100 (13.0)	100 (13.0)	100 (20.0)
<i>M. scrofulaceum</i>	1	100 (15.0)	100 (15.0)	100 (15.0)
<i>M. chelonae</i>	5	100 (32.0)	100 (32.0)	100 (23.8)

Limitations of the procedure

Mycobacteria with special growth requirements, such as *M. haemophilum* or *M. paratuberculosis* and other mycobactin-dependent mycobacteria, do not grow on solid media without supplementary nutrients. Furthermore, certain mycobacterial species, such as *M. genavense*, cannot be cultured on solid media⁹.

Quality control

Quality control tests are performed on each lot of Mycobacterium tube 1 and 2 at the time of manufacture. Should the user wish to perform his own quality control, the following organisms are recommended:

1. *M. tuberculosis* ATCC 25177
2. *M. intracellulare* ATCC 13950
3. *M. scrofulaceum* ATCC 19981
4. *M. kansasii* ATCC 12478
5. *M. fortuitum* ATCC 6841

Disposal

All inoculated tubes may contain highly contagious material and should be disposed of according to existing local regulations.

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1. Katila M-L, Mattila J. Enhanced isolation of MOTT on egg media of low pH. APMIS 1991;99:803–7.
2. Katila M-L, Mattila J, Brander E. Enhancement of growth of Mycobacterium malmøense by acidic pH and pyruvate. Eur J Microbiol Infect Dis 1989;8:998–1000.
3. Naumann L. Vergleichende Untersuchungen von flüssigen und festen Nährböden in der Tuberkulose-Diagnostik. Klinisches Labor 1992;38:325–9.
4. Portaels F, Pattyn SR. Growth of mycobacteria in relation to the pH of the medium. Ann Microbiol (Paris) 1982;133:212–21.
5. Stonebrink B. Tubercle bacilli and pyruvic acid. Proceed. Tub. Res. Council. 1957;44:67–74.
6. Shinnick TM, Good RC. Diagnostic mycobacteriology laboratory practices. Clin Infect Dis 1995;21:291–9.
7. Falkinham JO. Epidemiology of infection by nontuberculous mycobacteria. Clin Microbiol Rev 1996;9:177–215.
8. Lue YA, Brophy G, Visot B. New solid medium for the rapid recovery of mycobacteria. In: Proceedings of the 95th ASM General Meeting; 1995 May 21–25; Washington, DC, USA. 1995.
9. Nolte SF, Metchock B: Mycobacterium. Manual of clinical microbiology. 6th ed., 1995.

Explanation of symbols • Forklaring på symboler • Symbolien selitykset			
	English	Norsk	Suomi
	For <i>in vitro</i> diagnostic use	For <i>in vitro</i> diagnostikk	<i>In vitro</i> -diagnostiikkaan
	Catalogue number	Artikkel nummer	Tuotenumero
	Batch code	Batch nr.	Eränumero
	Use by	Brukes til	Käytettävä ennen
	Store at 2...8°C	Storage temperature	Oppbevarings-temperatur
	Manufacturer	Produsent	Valmistaja

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