Intended use

A dip slide culture method for diagnosing urinary tract infections by demonstrating microbes in urine.

Contents

Uricult Trio	Cat. No. 68197	Cat. No. 68300
Dip slides	10	300
Patient labels	10	300
Instructions for use	1	1

Principle

The Uricult Trio dip slide system is based on three agar media. One side of the plastic slide is covered with green CLED medium and the other with both reddish-brown MacConkey medium and a colourless E. coli medium.

The CLED medium is intended for determining the total bacterial count. The MacConkey medium is intended for detecting gram-negative organisms. This medium contains bile salts that prevent the growth of gram-positive organisms other than enterococci which may grow as minute colonies. The E. coli medium is intended especially for detecting gram-negative, β -glucuronidase-producing organisms which grow as colonies in various shades of brown or grey on this medium. *Escherichia coli* is the most common β -glucuronidase-producing organism encountered in urinary tract infections. Bile salts contained in the E. coli medium prevent the growth of gram-positive organisms on it. Certain yeasts do grow on the E. coli medium, however.

Typical formulation

CLED medium	MacConkey medium	E. coli medium
Peptone 10.0 g/l	Peptone 20.0 g/l	Peptone 12.0 g/l
Meat extract 3.0 g/l	Lactose 10.0 g/l	MgSO ₄ 0.1 g/l
Lactose 10.0 g/l	Neutral red 0.075 g/l	MnCl ₂ 0.01 g/l
L-Cystine 0.13 g/l	Bile salts 0.8 g/l	Fe(III) citrate 0.4 g/l
Bromthymol blue 0.03 g/l		Bile salt mixture 2.4 g/l
		8-hydroxyquinoline β -D glucuronide

Warnings and precautions

Uricult Trio is for *in vitro* diagnostic use only. Do not use the product beyond the expiry date marked on the box. Do not use the Uricult Trio if you detect discoloration or dehydration of the agar, separation of the growth media from the plastic slide or evidence of bacterial or fungal growth. Because any colonies growing on Uricult Trio are actual or potential pathogens, do not touch the growth.

Storage

Uricult Trio is stored at 7...25°C, protected from air and temperature fluctuations. Avoid draughts and storage near heat-generating appliances. Do not allow to freeze. The expiry date is marked on the box.

Urine sampling and the storage of samples

Ideally, urine for bacterial culture should remain in the bladder for four hours prior to sampling. Samples are usually obtained from clean-voided (mid-stream) urine. Samples for culture may also be obtained by catheterisation or suprapubic bladder aspiration.

The sample should be inoculated onto the Uricult Trio slide immediately after collection. The slide should then be returned into its protective tube and the cap closed tightly. If the urine sample needs to be stored prior to inoculation, it should be maintained refrigerated at 2...8°C no longer than 24 hours. Medications (especially antibiotics) may affect the result of the Uricult Trio test. Therefore, the test should not be performed until 48 hours after the final dose of medication.

Test procedure

- Unscrew the slide from the tube without touching the agar surfaces.
- Holding Uricult Trio by the cap, dip the slide into freshly voided mid-stream urine so that the agar surfaces are totally immersed. If the volume of urine is too small for this, the agar surfaces can be wetted by pouring urine on them, followed by tilting to ensure complete wetting.
- Allow excess urine to drain from the slide.
- Blot the last drops on absorbent paper.
- Screw the slide tightly back into the tube.
- Fill in the patient label and attach it to the tube.
- Place the tube upright in an incubator (36±2°C) for 16–24 hours. The tube may also be sent to a laboratory for incubation.
- To obtain a colony count (CFU/ml), remove the slide from the tube and compare the colony density with the model chart provided in the kit.

Note:

- Negative cultures may be incubated for additional 24 hours to ensure that slow-growing bacteria are detected.
- The inoculated slide may be incubated immediately or stored or transported to a laboratory for incubation and interpretation. Storage or transportation should not exceed 48 hours at 7...25°C, after which Uricult Trio should be incubated at 36±2°C for 16–24 hours. If the slide has been stored or transported for up to 48 hours, only the presence of growth and the colony count should be recorded from it; the colour reaction may be atypical.
- The inoculated slide may be incubated at room temperature for 1–3 days, after which positive cultures may be sent to a specialised laboratory for further investigation⁵. Negative cultures may be incubated for additional 24 hours to detect slow-growing bacteria⁶.

Interpretation of results

After incubation of the inoculated slide, the presence of bacteria is evidenced by colonies on the agar surface. Because a colony is the result of the multiplication of a single bacterial cell, the number of colonies indicates the concentration of colony-forming units (CFUs/ml) in the urine sample. The colony count should be determined from the originally green CLED medium by matching the colony density with the model chart it most closely resembles. It is important to compare the number of colonies, not their size.

The low electrolyte concentration of the CLED medium prevents spreading of Proteus strains. Bromthymol blue and lactose in the medium allow the detection of lactose-fermenting bacteria. Such lactose-positive strains grow as yellow colonies and turn the medium yellow, whereas lactose negative strains grow as translucent colonies with no colour change of the medium.

The originally brownish-red, selective MacConkey medium supports the growth of gram-negative bacteria, but even enterococci may grow as pinpoint colonies on the medium⁷. Lactose-positive bacteria grow as red and lactose-negative bacteria as translucent colonies on the medium. β -glucuronidase-producing organisms grow as colonies in various shades of brown or grey on the colourless E. coli medium^{8,9}. β -glucuronidase-negative, gram-negative bacteria grow as translucent colonies on this medium.

When the urinary bacterial content is high ($\geq 10^7$ CFU/ml), the agar surfaces may become totally covered by confluent growth. This can be misinterpreted as a negative result. Therefore, any surfaces that appear negative should be examined under a reflecting light; absence of reflection indicates confluent growth. A bright light also allows very small colonies to be detected.

A mixture of different bacterial strains on the Uricult Trio is most likely due to contamination of the urine sample.

Expected values

The following values are based on the ECLM-EUG European Urinalysis Guidelines (2000).

Method of sampling clinical status	Significant colony count (CFU/ml)
Midstream, bladder time < 4 hours, symptomatic patient	$\geq 10^3$
Midstream, bladder time > 4 hours	$\geq 10^{4-5}$
Catheter sample from man	$\geq 10^3$
Catheter sample from woman	$\geq 10^4$
Nonsymptomatic bacteriuria	$\geq 10^5$
Puncture sample	Any growth

Note: In some cases, bladder urine < 4 hours may express clinically significant colony counts below 10^3 CFU/ml.

Limitations of procedure

Uricult Trio is capable of detecting bacterial concentrations between 10^3 and 10^7 CFU/ml. The model chart allows the determination of colony counts to the nearest power of 10. When the chart is used according to instructions, colony counts show a 99 % correlation with the conventional pour plate method¹.

Performance characteristics**Uricult Trio • CLED medium**

Arnell, G.C. 1970: Detection of bacteriuria at room temperature. Lancet, January 17, pp 119–121 ¹ .		
Number of samples	140	Reference method: Pour plate (Nutrient agar)
Sensitivity	100 %	
Specificity	99 %	
PPV	98 %	
NPV	100 %	

Quality control

Quality control tests are performed on each lot of Uricult Trio dip slides at the time of manufacture. Should the user wish to perform his own quality control, the following procedure is recommended:

- Prepare a 10^5 – 10^8 bacteria/ml suspension of each of the following bacteria in sterile saline:
 - Staphylococcus aureus* ATCC 25923
 - Escherichia coli* ATCC 25922
 - Proteus mirabilis* ATCC 12453
- Use the suspensions to inoculate the Uricult Trio dip slides, using the normal procedure.
- Interpret the results after a 16–24 hour incubation as follows:


S. aureus ATCC 25923: Growth of colonies on the CLED medium only. Colonies ferment lactose, as indicated by the yellow colour of the medium and the shift towards yellow of the medium.

E. coli ATCC 25922: Growth of yellow colonies with a shift towards yellow of the CLED medium and growth of pink-red colonies on the MacConkey medium. Growth of brown colonies on the E. coli medium.


P. mirabilis ATCC 12453: Growth of translucent colonies on the CLED medium, with a possibility of bluish coloration of the medium. Growth of colourless colonies on the MacConkey medium. Growth of colourless colonies on the E. coli medium.

Disposal

Used Uricult Trio dip slides are best disposed of by burning, autoclaving or immersing in a disinfectant overnight, with adherence to local regulations.

Model Chart • Tableau de référence • Standardbildkarte • Tablas de referencia • Tabela de Referência • Tavola di riferimento • Πρότυπος πίνακας αναφοράς • Modelová tabulka • Referenčná tabela • Modelkort • Avlesningsmal • Tolkningsmall • Mallitaulu

Literature • Bibliographie • Literatur • Bibliografía • Referências Bibliográficas • Bibliografia • Βιβλιογραφία • Literatura • Referenser • Litteratur • Kirjallisuus

- McAllister TA, Arnell GC, Barr W, Kay P: Assessment of plain dipslide quantitation of bacteriuria. Nephron 11: 111–122, 1973.
- Kass EH: Bacteriuria and the diagnosis of infections of the urinary tract. Archives of Internal Medicine 100: 709–714, 1957.
- Mackey JP, Sandys GH: Laboratory diagnosis of infections of the urinary tract in general practice by means of a dip-inoculum transport medium. British Medical Journal 2: 1286–1288, 1965.
- NCCLS Publication M22-A: Quality Assurance Standards for Commercially Prepared Microbiological Culture Media; Approved standard. Volume 10: 14, 1990.
- Ekman et al.: Streptococci- ja Uricult-levyjen käyttö päivystysaikana. Aesculapius 11–12, 1985.
- Arnell GC: Detection of bacteriuria at room temperature. Lancet, January 17: 119–121, 1970.
- Granato PA: Evaluation of a dip slide device for enumeration of bacteria in urine. Laboratory Medicine Vol. 11, No 4: 246–250, 1980.
- Dalet F, Segovia T: Evaluation of a New Agar in Uricult Trio for Rapid Detection of *Escherichia coli* in Urine. J. Clin. Microbiol. 1395–1398, 1995.
- Larinkari U, Rautio M: Evaluation of New Dipslide with a Selective Medium for the rapid detection of Beta-Glucuronidase-Positive *Escherichia coli*. Eur. J. Clin. Microbiol. Infect. Dis. Vol. 14: 606–608, 1995.

Explanation of symbols • Explication des symboles • Erläuterung der Symbole • Explicación de los símbolos • Explicação de símbolos • Spiegazione dei simboli • Επεξήγηση των συμβόλων • Vysvětlivky použitých symbolů • Pojasnila simbolov • Forklaring af symboler • Förklaring av symboler • Symbolien selitykset

IVD	REF	LOT	Store at 7...25°C		
For <i>in vitro</i> diagnostic use Diagnostic <i>in vitro</i> In Vitro Diagnostikum Para uso diagnóstico <i>in vitro</i> Para uso em diagnóstico "in Vitro" Per uso diagnostico <i>in vitro</i> Για διαγνωστική χρήση μόνο Za <i>in vitro</i> diagnostiko Kun til <i>in vitro</i> diagnostikk brug For <i>in vitro</i> diagnostikk För <i>in vitro</i> diagnostik In vitro -diagnostiikkaan	Catalogue number Référéncé Bestell Nr. Número de catálogo Número de catalogo (no) Codice Αριθμός καταλόγου Katalogovné číslo Kataloška številka Varennummer Artikkel nummer Artikelnummer Tuotenumero	Batch code Número de lot Lot. Nr. Número de lote Número de lote Numero di lotto Αριθμός παρτίδας Císlo sarže Serijaska številka Lot nr. Batch nr. Batch nr. Eräkoodi	Storage temperature Conservation Lagertemperatur Conservación Conservação Conservazione Θερμοκρασία φύλαξης Teplotní rozmezí Temperatura sхранjenja Opbevaringstemperatur Oppbevaringstemperatur Förvaringstemperatur Säilytyslämpötila	Use by a utiliser jusqu'à Verwendbar bis Fecha de caducidad Utilizaz até Utilizzare entro/Scadenza Χρησιμοποιήστε έως Použitelné do Uporabno do Anvendes inden Brukes til Används före Käytettävä ennen	Manufacturer Fabricant Hersteller Fabricante Produtor Produttore Κατασκευάζεται από Výrobce Pouzavajalec Producent Tilværkare Valmistaja

Cap • Bouchon • Deckel • Tapón • Tampa
Tappo • Καπάκι • Víčko • Pokrov • Låg
Lock • Kanski

Slide • Lame • Träger • Lámina • Placa
Plastra • Flávka • Destička • Testna ploščica
Slida • Levy

Tube • Röhrchen • Tubo • Cilindro
Συλινδρος • Nädobka • Epruvetka
Rör • Rör • Putki

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