Intended use

A selective enrichment broth for primary isolation and cultivation of *Trichomonas vaginalis* and yeasts from urogenital samples. *T. vaginalis* and yeasts are common causative agents of vaginitis¹. *T. vaginalis* is also an inhabitant of the male urogenital system².

Principle of the method

The 4 ml trypticase-based nutrient broth in the Vagicult tube is a modification of Diamonds medium^{3,4}. It contains sheep serum, with penicillin and streptomycin serving as selective agents to suppress the growth of contaminating bacteria.

Contents

Vagicult	Cat. No. 67693
Vagicult enrichment broth tubes	20 pcs
Patient labels	20 pcs
Instructions for use	1 pc

Material required but not included

Laboratory equipment for microbial culture and microscopy with x40 phase-contrast and/or differential interference contrast optics⁵.

Typical formulation

Trypticase peptone	22.0 g/l	Dihydrostreptomycin	
Yeast extract	11.0 g/l	sulphate	1.0 g/l
Maltose	5.5 g/l	Penicillin G	0.6 g/l
Cysteine	1.1 g/l	Sheep serum	100 ml/l
Ascorbic acid	0.2 g/l		

Warnings and precautions

- Vagicult is for **in vitro diagnostic use** only.
- Do not use the product beyond the expiry date.
- Do not use the Vagicult tube if you detect evidence of bacterial or fungal growth or if the tube is turbid.
- Because any growth in Vagicult tubes is potentially pathogenic, do not touch it.

Storage

Vagicult tubes are stored refrigerated (2...8°C). The expiry date is marked on the box.

During storage, some sediment may accumulate on the bottom of the tube, appearing as fibrel-like material or small floccules upon shaking. This will not impair the performance of the tube.

Sample collections

Before inoculation, warm the Vagicult tube to 25...37°C.

1. Inoculate the Vagicult tube with the patient sample.

Women:

- a. Take a sample with a swab from vaginal discharge or the posterior fornix.
- b. Transfer the sample immediately to a Vagicult tube by rotating the swab vigorously in the enrichment broth.
- c. Do not leave the swab in the broth.

Men:

- a. Inoculate the Vagicult tube with a drop of semen, prostatic fluid, urethral swab or concentrated urine.
 - b. Rotate the tube vigorously.
2. Fill in the patient label and attach it to the tube.
 3. Place the tube upright in an incubator at 36±2°C.
 4. Incubate the tube up to 72 h.

Reading the result

Both *T. vaginalis* and yeasts make the broth turbid during growth. A sample is taken daily from the growth on a glass slide, and three visual fields are examined under a microscope at x40 magnification.

T. vaginalis can be seen as a free-swimming, pear-shaped parasite. Yeast cells are round or oval, often showing budding.

The total number of *T. vaginalis* parasites or yeast cells in the three visual fields is interpreted as follows:

1–5 pcs:	+
6–30 pcs:	++
over 30 pcs:	+++

After a 24-h incubation, a yeast result of ++ is significant. Additional identification for *Candida* can be performed using the Mackenzie serum tube test or at a specialised laboratory.

After a 24-h incubation, a *T. vaginalis* result of + is significant.

Expected values

T. vaginalis inhabits the urogenital systems of both men and women. The parasite is considered a pathogen. The infection in men is typically asymptomatic, whereas 25–50 % of infected women exhibit symptoms, such as dysuria, vaginal itching or burning. In severe infections, vaginal discharge is encountered².

Symptomatic vaginal candidosis has been experienced by 40–75 % of sexually active women².

Performance characteristics

Culturing is the most sensitive method of diagnosing trichomoniasis. Along with culture, wet smears and/or stained smears should be examined microscopically^{2,5,6}.

T. vaginalis

McMillan A. Laboratory diagnostic methods and cryopreservation of Trichomonads. In: B.M. Honigberg (ed.) Trichomonads parasitic in Man, Springer-Verlag, New York, 1990⁶. (Summary reference of five separate studies using different *T. vaginalis* culture media.)

Comparison of percentages of <i>T. vaginalis</i> infections diagnosed by wet smear and/or culture			
Total no. of samples	Wet smear	Culture	Combination
100	74	83	91
107	48	86	100
175	86	78	100
198	88	79	98
54	52	89	93

Candida

Carlson P, Richardson M, Paavonen J: Evaluation of the Oricult N dipslide for laboratory diagnosis of vaginal candidiasis. J. Clin. Micro. 38 (3): 1063–1065, 2000.

Proportions of vaginal swabs positive for vaginitis by microscopy and culture in vaginitis patients and control subjects

Group of samples	Total no.	Microscopy	Vagicult
Patients	124	14 (11 %)	21 (17 %)
Controls	50	2 (4 %)	11 (22 %)

Limitations of the procedure

A broth culture method does not differentiate between *Candida* patients with disease and asymptomatic carriers⁷.

Because the culture method may take up to three days, and the patient sample may contain nonviable organisms, it is important that wet smears are also examined⁵.

Quality control

Quality control tests are performed on each lot of Vagicult at the time of manufacture. Should the user wish to perform his own quality control, the following organisms are recommended:

Candida albicans

1. Warm the tubes at 35...37°C for at least 20 min.
2. Inoculate a tube with a sterile loop carrying one small colony of *C. albicans* ATCC 8014.
3. Inspect a wet mount from the broth under a microscope (x40 magnification) and count the cells in 3 visual fields.
4. Incubate the tube at 35...37°C for 48 h.

T. vaginalis

5. Inoculate a tube with 5–8 drops from a 48-h-old or older stock culture of living *T. vaginalis* ATCC 30001 cells.
6. Incubate the tube at 35...37°C for 48 h.
7. Inspect a wet mount from the broth under a microscope (x40 magnification) and count the cells in 3 visual fields.
8. After incubation, repeat steps 3 and 7 to confirm an increase in the number of cells in both tubes.

Disposal

Inoculated Vagicult tubes with detectable growth may contain highly contagious material. Inoculated tubes should be disposed of according to existing local regulations.

Vagicult is a registered trademark of Orion Diagnostica Oy.

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